



Default Layout

Table EXPERIMENT_RECORDS			
* Pk	record_id	integer	
*	LOT_NO	integer	
*	project_title	varchar(256)	
*	formulation_date	date	
*	preparation_date	date	
*	prepared_by	varchar(256)	
*	quantity	double(10,2)	Quantity in grams
*	total_percentage_w	double(10,2)	
*	total_AR	double(10,2)	
*	total_AD	double(10,2)	
	notes	varchar(600)	
	preparation_reason	varchar(600)	
	observations	varchar(600)	
	obs_description	varchar(256)	Single word description of experiment (e.g transparent)
	rating	varchar(256)	rating of experiment (good/bad)
*	date_created	datetime	
	date_updated	datetime	
Indexes			
Pk	Pk_EXPERIMENT_RECORD	record_id	
	S_record_id		

Table HYDROGEN_PEROXIDE_DATA			
list of hydrogen peroxide and hydrogen peroxide stability results for the experiment records			
* Pk	hp_id	integer	
*	experiment_record_id	integer	
*	hp_type	integer	records the type of HP data recorded (hydrogen peroxide or hydrogen peroxide stability)
*	trial_no	integer	
*	N	double(10,4)	
	M	double(10,4)	
	vol_change	double(10,4)	change in volume in ml
*	H2O2	double(10,4)	
	PH	double(10,4)	
	accepted_range	varchar(50)	
	date	date	
	initials	varchar(50)	
	conditions	varchar(256)	
	SG	varchar(256)	
Indexes			
Pk	Pk_HYDROGEN_PEROXIDE_DATA_hp_id	hp_id	
Foreign Keys			
	Fk_HYDROGEN_PEROXIDE_DATA_record ( experiment_record_id ) ref EXPERIMENT_RECORDS ( record_id )		

Table RAW_MATERIALS			
list of all raw materials			
* Pk	raw_material_id	integer	
*	experiment_record_id	integer	foreign key which links the raw materials to the record_id from EXPERIMENT_RECORDS table
*	percentage_w	double(10,2)	
	raw_material_lot	integer	
*	AR	double(10,2)	

Table RAW\_MATERIALS

*	AD	double(10,2)	
	notes	varchar(600)	
Indexes			
Pk	Pk_RAW_MATERIALS_raw_material_id	raw_material_id	
Foreign Keys			
	Fk_RAW_MATERIALS_id ( experiment_record_id ) ref EXPERIMENT_RECORDS ( record_id )		